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
BONE MARROW BIOPSY NEEDLE**INSTRUCTIONS FOR USE:****ENG****Indications For Use:** Intended for the purpose of harvesting bone marrow specimens.**Possible Side Effects:** Pain, fever, bleeding, hematuria and infection.**Caution:** Ethylen oxide sterilized. Non-pyrogenic. Check the package. if package damaged or opened do not use the products. Check the expiry date and the gauge. Store in dry place, between 5 – 30 °C temperature, protect from sun light. Possible allergic reactions should be considered. It must be used only in hospitals by physicians who familiar with the device and procedures. For single Patient Use Only. Do not attempt to clean or resterilize this product. Do not use the device under MRI. After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.**Please double check for device function and possible damage for correct implantation and do not use if package is opened or damaged.****ASPIRATION PROCEDURE:**

- 1.Place the patient in a right or left lateral position, with comfortably flexed and the top knee drawn toward the chest.
- 2.Locate the posterior superior iliac spine and mark.
- 3.Using the sterile technique, prepare the skin with antiseptic and drape.
- 4.Infiltrate the marked area with local anaesthetic, especially the periosteum.
- 5.Make a skin incision with a scalpel blade over the marked area.
- 6.Hold the needle with proximal end in palm and the index finger against the shaft near the tip. This position stabilizes the needle and allows better control.
- 7.Introduce the needle through the incision pointing toward the anterior superior iliac spine and bring it into contact with the posterior iliac spine.
- 8.Using gentle, but firm pressure, advance the needle. Rotate the needle in an alternating clockwise - counterclockwise motion. The entrance into the marrow cavity is generally detected by a decreased resistance.
- 9.Remove the stylet by pulling the upper section of the handle straight out.
- 10.Attach a syringe with a luer taper to the hub of the biopsy needle with a firm, push, twist motion.
- 11.Apply suction by withdrawing syringe plunger. Remove the aspirated specimen collected from the syringe.

BIOPSY PROCEDURE:

- Take the biopsy after the aspiration using the same skin incision, but choosing a different location to obtain the biopsy sample.
- Repeat steps 1 through 9 of the procedure and continue as follows:
10. Optional - If desired, cannula can be capped with luer cap provided.
 11. Slowly and gently advance the needle a millimeter at a time, with clockwise - counterclockwise motion until adequate marrow is obtained.
 12. Remove the needle by pulling it straight out. Insert the probe to the needle through the tip adapter in order to push the specimen out.

Note: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

STERILE EO

Single Use Storage limits Expiry Date Ethylen oxide sterilized

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GEOTEK MEDİKAL VE SAĞLIK HİZ.LTD.STİ