

BONE MARROW BIOPSY NEEDLE

KEMİK İLİĞİ BIYOPSİ İĞNESİ

- ↳ Atraumatic ultra sharp needle tip design
- ↳ Cutting Cannula helps to easier and faster sampling
- ↳ Light weight and ergonomic handle
- ↳ Easy to use with excellent comfort



Ref. No	Diameter (G)	Length (cm)
GBB0810	8	10
GBB0815	8	15
GBB0910	9	10
GBB0915	9	15
GBB1110	11	10
GBB1115	11	15
GBB1310	13	10
GBB1315	13	15

- ↳ Atravmatik ultra keskin uç dizayn
- ↳ Hafif ve ergonomik sap
- ↳ Kolay ve hızlı numune alımı sağlayan kesme kanülü
- ↳ Konforlu ve kolay uygulama

w / Cutting Canula

Ref. No	Diameter (G)	Length (cm)
GBB0810-AL	8	10
GBB1110-AL	11	10

