

# DOUBLE J URETERAL STENT SET


## DOUBLE J ÜRETERAL STENT SET

- ↳ Soft polyurethane material
- ↳ Higher kink resistance
- ↳ High radiopacity
- ↳ Depth marker
- ↳ Long pusher option between 40 cm - 100 cm
- ↳ Suture on distal tip



| Diameter (F) | Length (cm) | Ref. No: (One End Open) | Ref. No: (Both End Open) |
|--------------|-------------|-------------------------|--------------------------|
| 3            | 12          | GDJ312-OEC              | GDJ312-BEO               |
|              | 14          | GDJ314-OEC              | GDJ314-BEO               |
|              | 16          | GDJ316-OEC              | GDJ316-BEO               |
|              | 20          | GDJ320-OEC              | GDJ320-BEO               |
| 4            | 12          | GDJ412-OEC              | GDJ412-BEO               |
|              | 14          | GDJ414-OEC              | GDJ414-BEO               |
|              | 16          | GDJ416-OEC              | GDJ416-BEO               |
|              | 20          | GDJ420-OEC              | GDJ420-BEO               |
|              | 22          | GDJ422-OEC              | GDJ422-BEO               |
|              | 24          | GDJ424-OEC              | GDJ424-BEO               |
| 4.7          | 26          | GDJ426-OEC              | GDJ426-BEO               |
|              | 12          | GDJ4712-OEC             | GDJ4712-BEO              |
|              | 16          | GDJ4716-OEC             | GDJ4716-BEO              |
|              | 20          | GDJ4720-OEC             | GDJ4720-BEO              |
|              | 22          | GDJ4722-OEC             | GDJ4722-BEO              |
|              | 24          | GDJ4724-OEC             | GDJ4724-BEO              |
| Adjustable   | 26          | GDJ4726-OEC             | GDJ4726-BEO              |
|              | 28          | GDJ4728-OEC             | GDJ4728-BEO              |
|              | 22/30       | VDJ47-OEC               | VDJ47-BEO                |
|              | 4.8         | 12                      | GDJ4812-OEC              |
| 16           |             | GDJ4816-OEC             | GDJ4816-BEO              |
| 20           |             | GDJ4820-OEC             | GDJ4820-BEO              |
| 22           |             | GDJ4822-OEC             | GDJ4822-BEO              |
| 24           |             | GDJ4824-OEC             | GDJ4824-BEO              |
| 26           |             | GDJ4826-OEC             | GDJ4826-BEO              |
| Adjustable   | 28          | GDJ4828-OEC             | GDJ4828-BEO              |
|              | 22/30       | VDJ48-OEC               | VDJ48-BEO                |
| 5            | 16          | GDJ516-OEC              | GDJ516-BEO               |
|              | 20          | GDJ520-OEC              | GDJ520-BEO               |
|              | 22          | GDJ522-OEC              | GDJ522-BEO               |
|              | 24          | GDJ524-OEC              | GDJ524-BEO               |
|              | 26          | GDJ526-OEC              | GDJ526-BEO               |
|              | 28          | GDJ528-OEC              | GDJ528-BEO               |
| 6            | 20          | GDJ620-OEC              | GDJ620-BEO               |
|              | 22          | GDJ622-OEC              | GDJ622-BEO               |
|              | 24          | GDJ624-OEC              | GDJ624-BEO               |
|              | 26          | GDJ626-OEC              | GDJ626-BEO               |
|              | 28          | GDJ628-OEC              | GDJ628-BEO               |
|              | 30          | GDJ630-OEC              | GDJ630-BEO               |
| Adjustable   | 22/30       | VDJ06-OEC               | VDJ06-BEO                |
|              | 7           | 22                      | GDJ722-OEC               |
| 24           |             | GDJ724-OEC              | GDJ724-BEO               |
| 26           |             | GDJ726-OEC              | GDJ726-BEO               |
| 28           |             | GDJ728-OEC              | GDJ728-BEO               |
| 8            | 30          | GDJ730-OEC              | GDJ730-BEO               |
|              | 22          | GDJ822-OEC              | GDJ822-BEO               |
|              | 24          | GDJ824-OEC              | GDJ824-BEO               |
|              | 26          | GDJ826-OEC              | GDJ826-BEO               |
| Adjustable   | 28          | GDJ828-OEC              | GDJ828-BEO               |
|              | 30          | GDJ830-OEC              | GDJ830-BEO               |

- ↳ Yumuşak poliüretan materyal
- ↳ Yüksek bükülme direnci
- ↳ Yüksek radyopasite
- ↳ Derinlik ölççeği
- ↳ 40-100 cm arasında uzun itici seçeneği
- ↳ Radyopak itici seçeneği
- ↳ Alt uçta suture



### DJ + PUSHER ONLY

Please Use Ref. No as GDJPPX-OEC or BEO

### HYDROPHILIC DJ STENT SET

Please Use Ref. No as GDJXX-HOEC or HBEO

- ↳ Set includes DJ stent + Pusher + Guide wire + Two Clamp