

# KATRACORE® AUTOMATIC BIOPSY NEEDLE

## KATRACORE® OTOMATİK BİYOPSİ İĞNESİ

- ↳ Easy and one hand operated with two stage arming
- ↳ Echogenic and atraumatic needle tip design
- ↳ 22mm advancement with 20mm notch
- ↳ Ergonomic and light weight design
- ↳ Depth marked canula
- ↳ Coaxial needles are optionally available



Ref. No	Diameter (G)	Length (cm)	Coaxial Needle Length (cm)	With Coaxial Ref. No
GKR1410	14	10	5	GKR1410-X
GKR1415	14	15	10	GKR1415-X
GKR1420	14	20	15	GKR1420-X
GKR1610	16	10	5	GKR1610-X
GKR1615	16	15	10	GKR1615-X
GKR1620	16	20	15	GKR1620-X
GKR1810	18	10	5	GKR1810-X
GKR1815	18	15	10	GKR1815-X
GKR1820	18	20	15	GKR1820-X
GKR1825	18	25		
GKR1830	18	30		
GKR2010	20	10	5	GKR2010-X
GKR2015	20	15	10	GKR2015-X
GKR2020	20	20	15	GKR2020-X

- ↳ İki aşamalı kurma sistemi sayesinde kolay ve tek elle kurulma
- ↳ Ekojenik ve atravmatik uç yapısı
- ↳ 22mm doku içine gömülme mesafesi ve 20mm örnek alma haznesi
- ↳ Ergonomik ve hafif sap dizaynı
- ↳ Derinlik işaretli kanül
- ↳ Koaksiyel iğneler isteğe bağlı sağlanabilir